



Product Datasheet

PTCH1 antibody (orb125007)


Description	Goat polyclonal antibody to PTCH1
Species/Host	Goat
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Target	PTCH (C Terminus)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	IHC: Paraffin embedded Human Brain (Coretx). Recommended concentration: 8µg/ml. Notes: Immunofluorescence: Strong expression of the protein seen in the Golgi of HeLa and NIH3T3 cells. Recommended concentration: 10µg/ml.
Protein Sequence	ELQDVECEERPR
Clonality	Polyclonal
MW	161; 161; 154; 144
Entrez	5727
Dilution Range	ELISA: 1:16000, IHC-P: 8 µg/ml
Expiration Date	12 months from date of receipt.



orb125007
(3.75µg/ml) staining
of paraf...



Immunofluorescence
analysis of HeLa
cell...



orb125007
(3.8µg/ml) staining
of paraff...