

Product Datasheet

ROBO3 antibody (orb125008)

Description	Goat polyclonal antibody to ROBO3
Species/Host	Goat
Reactivity	Canine, Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Target	ROBO3 (internal)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Notes: Flow Cytometry: Flow cytometric analysis of U2OS cells. Recommended concentration: 10ug/ml. This antibody was successfully used on rheumatoid arthritis synovial fibroblasts (RASf): Denk et al (2010) PMID: 20298552.
Protein Sequence	RKDGARLKEEEGR
Clonality	Polyclonal
MW	148
Entrez	64221
Dilution Range	ELISA: 1:64000, IHC-P: 10ug/ml
Expiration Date	12 months from date of receipt.