


Product Datasheet

SFRP1 antibody (orb125011)

Description	Goat Polyclonal to SFRP1. NP_003003.3, PMID: 17699851 and GenElD: 6422.
Species/Host	Goat
Reactivity	Bovine, Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Target	SFRP1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	IHC: In paraffin embedded Human Kidney shows vesicular staining in epithelial cells of PCT. Recommended concentration, 3-6µg/ml.
Protein Sequence	KNHECPTFQSVFK
Clonality	Polyclonal
MW	35.4
Entrez	6422
Dilution Range	ELISA: 1:1000, IHC-P: 3-6µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human Ki...