


Product Datasheet

PCSK9 antibody (orb125015)

Description	Goat Polyclonal to PCSK9. NP_777596.2, PMID: 18250299 and GeneID: 255738.
Species/Host	Goat
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, WB
Target	PCSK9 (aa214-228)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 70-75kDa band observed in HepG2 cell lysates and approx 75kDa in HeLa cell lysates (calculated MW of 74.3kDa according to NP_777596.2).. Recommended concentration: 1-2g/ml. Primary incubation 1 hour at room temperature.. Notes: Flow Cytometry: Flow cytometric analysis of HeLa cells. Recommended concentration: 10ug/ml.
Protein Sequence	TRFHRQASKCDSHGT
Clonality	Polyclonal
MW	74.3
Entrez	255738
Dilution Range	ELISA: 1:16000, WB: 1-2g/ml, IHC-P: 10ug/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293 lysate u...