

Product Datasheet

SLC11A2 Antibody (orb1250182)

Description	SLC11A2 Antibody
Species/Host	Rabbit
Reactivity	Human, Primate
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Synthetic 15 amino acid peptide from C-Terminus of human SLC11A2.
Target	SLC11A2
Preservatives	PBS, 0.1% sodium azide.
Concentration	1 mg/ml
Storage	Antibody may be stored short term at 4°C, long term aliquot and store at -20°C, avoid freeze thaw cycles.
Note	For research use only
Application notes	SLC11A2 Antibody can be used for Immunohistochemistry in paraffin embedded tissues starting at 15 µg/mL.
Clonality	Polyclonal
Uniprot ID	P49281
NCBI	P49281
Dilution Range	SLC11A2 Antibody can be used for Immunohistochemistry in paraffin embedded tissues starting at 15 µg/mL.
Expiration Date	12 months from date of receipt.