


Product Datasheet

POU2F1 antibody (orb125019)

Description	Goat polyclonal antibody to POU2F1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	POU2F1 / OCT1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 100kDa band observed in nuclear lysates of cell line HeLa (calculated MW of 76.5kDa according to NP_002688.2). This molecular weight is routinely observed by other sources. Recommended concentration: 0.5-2 µg/ml.
Protein Sequence	NVQSKSNEESGDSQQ
Clonality	Polyclonal
MW	76.5
Entrez	5451
Dilution Range	ELISA: 1:4000, WB: 0.5-2 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa lysate using POU2F1 antibody.