


Product Datasheet

POU4F3 antibody (orb125020)

Description	Goat polyclonal antibody to POU4F3
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	POU4F3
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 40kDa band observed in Human and Mouse Heart lysates (calculated MW of 37.1kDa according to Human NP_002691.1). Recommended concentration: 0.2-0.6µg/ml.
Protein Sequence	EAAYREKNSKPE
Clonality	Polyclonal
MW	37.1
Entrez	5459
Dilution Range	ELISA: 1:8000, WB: 0.2-0.6 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Heart lys...