

Product Datasheet

CD79B Antibody (orb1250218)

Description

CD79B Antibody

Species/Host	Rat
Reactivity	Canine, Human, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC, IP, WB
Immunogen	Peptide containing 20 amino acids residues from mouse CD79b.
Target	CD79B
Preservatives	Tris buffered saline, pH 8.0, 0.09% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/ml
Storage	CD79B Monoclonal antibody may be stored short term at 4°C, long term aliquot and store at -20°C, avoid freeze thaw cycles.
Note	For research use only
Application notes	CD79B Monoclonal Antibody can be used for flow cytometry starting at 1:50 - 1:100 and immunohistochemistry (tested on paraffin embedded tissues) starting at 2 µg/mL.
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P40259
NCBI	P40259
Dilution Range	CD79B Monoclonal Antibody can be used for flow cytometry starting at 1:50 - 1:100 and immunohistochemistry (tested on paraffin embedded tissues) starting at 2 µg/mL.
Expiration Date	12 months from date of receipt.