

Product Datasheet

NEFM Antibody (orb1250231)

Description

NEFM Antibody

Species/Host	Mouse
Reactivity	Frog, Guinea pig, Hamster, Human, Monkey, Rabbit
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	Neurofilament preparation of calf brain tissue
Target	NEFM
Preservatives	PBS containing 0.09% sodium azide
Form/Appearance	Liquid
Concentration	1 mg/ml
Storage	NEFM Monoclonal antibody may be stored short term at 4°C, long term aliquot and store at -20°C, avoid freeze thaw cycles. Store undiluted.
Note	For research use only
Application notes	NEFM Monoclonal Antibody can be used for flow cytometry starting at 1:50 - 1:100, immunohistochemistry (tested on frozen tissues starting at 1:50 - 1:100, immunohistochemistry (tested on parafin embedded tissues) starting at 10 µg/mL and Western Blot starting at 1:100 - 1:500.
Clonality	Monoclonal
Uniprot ID	P07197
NCBI	P07197
Dilution Range	NEFM Monoclonal Antibody can be used for flow cytometry starting at 1:50 - 1:100, immunohistochemistry (tested on frozen tissues starting at 1:50 - 1:100, immunohistochemistry (tested on parafin embedded tissues) starting at 10 µg/mL and Western Blot starting at 1:100 - 1:500.
Expiration Date	12 months from date of receipt.