


Product Datasheet

CCR5 Antibody (orb1250256)

Description	CCR5 Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC, IP, WB
Immunogen	Human CCR5
Target	CCR5
Preservatives	Aqueous buffer, 0.09% sodium azide. Sourced from tissue culture supernatant.
Form/Appearance	Liquid
Concentration	0.5 mg/ml
Storage	Store antibody at 4 °C.
Note	For research use only
Application notes	CCR5 antibody can be used in Western Blot starting at 1 - 3 µg/mL, and immunohistochemistry starting at 20 µg/mL.
Isotype	IgG1,k
Clonality	Monoclonal
Uniprot ID	P51681
NCBI	P51681
Dilution Range	CCR5 antibody can be used in Western Blot starting at 1 - 3 µg/mL, and immunohistochemistry starting at 20 µg/mL.
Expiration Date	12 months from date of receipt.



Immunohistochemistry
of human tonsil tis...