



Product Datasheet

GRM2 antibody (orb125026)

Description	Goat polyclonal antibody to GRM2
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	GRM2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	IHC: In paraffin embedded Human Cerebral Cortex shows membranous/axonal staining. Recommended concentration, 3-5µg/ml.
Protein Sequence	KETAPERREV
Clonality	Polyclonal
MW	95.6; 27.5
Entrez	2912
Dilution Range	ELISA: 1:16000, WB: 1-3 µg/ml, IHC-P: 3-5µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis
of human
Hippocamp...



Immunohistochemical
analysis of formalin...