


Product Datasheet

APOD Antibody (orb1250299)

Description	APOD Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Apolipoprotein D antibody was raised against a synthetic peptide of Apolipoprotein D (Mouse).
Target	APOD
Preservatives	PBS, 0.01% sodium azide, 50% glycerol.
Form/Appearance	Liquid
Concentration	0.5 mg/ml
Storage	Apolipoprotein D antibody can be stored short term 4 °C. For long term storage aliquot and store at -20 °C. As with all antibodies avoid freeze/thaw cycles.
Note	For research use only
Application notes	Apolimmunoprecipitationoprotein D antibody can be used in ELISA, Western Blot, and immunohistochemistry starting at 5 µg/mL.
Clonality	Polyclonal
Uniprot ID	P05090
NCBI	P05090
Dilution Range	Apolimmunoprecipitationoprotein D antibody can be used in ELISA, Western Blot, and immunohistochemistry starting at 5 µg/mL.
Expiration Date	12 months from date of receipt.



Immunohistochemistry staining of Apolipo...