

Product Datasheet

HDAC4 Antibody (orb1250329)

Description	HDAC4 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	Histone Deacetylase 4 antibody was raised against amino acids 194 - 209 of Histone Deacetylase 4 (Human).
Target	HDAC4
Preservatives	PBS, 0.05% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/ml
Storage	Histone Deacetylase 4 antibody can be stored short term 4 °C. For long term storage aliquot and store at -20 °C. As with all antibodies avoid freeze/thaw cycles.
Note	For research use only
Application notes	Histone Deacetylase 4 antibody can be used in immunohistochemistry starting at 10 µg/mL.
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P56524
NCBI	P56524
Dilution Range	Histone Deacetylase 4 antibody can be used in immunohistochemistry starting at 10 µg/mL.
Expiration Date	12 months from date of receipt.