


Product Datasheet

EPHA2 Antibody (orb1250387)

Description	EPHA2 Antibody
Species/Host	Rabbit
Reactivity	Bovine, Canine, Equine, Human, Monkey, Mouse, Primate, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A portion of amino acids 450-500 of human EphA2.
Target	EPHA2
Preservatives	PBS, 0.2% gelatin, 0.05 % sodium azide.
Form/Appearance	Liquid
Concentration	0.5 mg/ml
Storage	EPHA2 antibody can be stored short term 4 °C. For long term storage aliquot and store at -20 °C. As with all antibodies avoid freeze/thaw cycles.
Note	For research use only
Application notes	EPHA2 antibody can be used in ELISA, Western Blot, immunohistochemistry starting at 5 µg/mL, and immunofluorescence starting at 10 µg/mL.
Clonality	Polyclonal
Uniprot ID	P29317
NCBI	P29317
Dilution Range	EPHA2 antibody can be used in ELISA, Western Blot, immunohistochemistry starting at 5 µg/mL, and immunofluorescence starting at 10 µg/mL.
Expiration Date	12 months from date of receipt.



Immunohistochemistry staining of EPHA2 i...