


Product Datasheet

SYK Antibody (orb1250408)

Description	SYK Antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	SYK monoclonal antibody was raised against amino acids 5 - 360 of SYK (Human).
Target	SYK
Preservatives	Phosphate-buffered solution, pH 7.2, 0.09% sodium azide.
Form/Appearance	Liquid
Concentration	0.5 mg/ml
Storage	Store SYK antibody at 4 °C, Avoid Freezing
Note	For research use only
Application notes	SYK antibody can be used in ELISA, and immunohistochemistry starting at 10 µg/mL.
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P43405
NCBI	P43405
Dilution Range	SYK antibody can be used in ELISA, and immunohistochemistry starting at 10 µg/mL.
Expiration Date	12 months from date of receipt.



Immunohistochemistry staining of SYK in ...