


Product Datasheet

BACE1 Antibody (orb1250439)

Description	BACE1 Antibody
Species/Host	Rabbit
Reactivity	Bat, Bovine, Canine, Equine, Guinea pig, Hamster, Human, Monkey, Mouse, Porcine, Primate, Rabbit, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	BACE1 antibody was raised against a peptide located near the internal domain of BACE1 (Human).
Target	BACE1
Preservatives	PBS, 0.1% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/ml
Storage	BACE1 antibody should be stored long term (months) at -80 °C and short term (days) at 4 °C. As with all antibodies avoid freeze/thaw cycles.
Note	For research use only
Application notes	BACE1 antibody can be used in immunohistochemistry starting at 1:50.
Clonality	Polyclonal
Uniprot ID	P56817
NCBI	P56817
Dilution Range	BACE1 antibody can be used in immunohistochemistry starting at 1:50.
Expiration Date	12 months from date of receipt.



Immunohistochemistry staining of BACE1 i...