


Product Datasheet

MST1R Antibody (orb1250457)

Description	MST1R Antibody
Species/Host	Rabbit
Reactivity	Human, Primate
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	MST1R antibody was raised against a peptide located near the internal domain of MST1R (Human).
Target	MST1R
Preservatives	PBS, 0.1% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/ml
Storage	MST1R antibody should be stored long term (months) at -80 °C and short term (days) at 4 °C. As with all antibodies avoid freeze/thaw cycles.
Note	For research use only
Application notes	MST1R antibody can be used in immunohistochemistry starting at 20 µg/mL.
Clonality	Polyclonal
Uniprot ID	Q04912
NCBI	Q04912
Dilution Range	MST1R antibody can be used in immunohistochemistry starting at 20 µg/mL.
Expiration Date	12 months from date of receipt.



Immunohistochemistry staining of MST1R i...