


Product Datasheet

GPR88 Antibody (orb1250472)

Description	GPR88 Antibody
Species/Host	Rabbit
Reactivity	Bovine, Gorilla, Hamster, Human, Monkey, Mouse, Rabbit, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P
Immunogen	GPR88 antibody was raised against a peptide located in the 1st cytoplasmic domain of GPR88 (Human).
Target	GPR88
Preservatives	PBS, 0.1% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/ml
Storage	GPR88 antibody should be stored long term (months) at -80 °C and short term (days) at 4 °C. As with all antibodies avoid freeze/thaw cycles.
Note	For research use only
Application notes	ELISA, IHC-P (10 - 20 ug/ml)
Clonality	Polyclonal
Uniprot ID	Q9GZN0
NCBI	Q9GZN0
Dilution Range	ELISA, IHC-P (10 - 20 ug/ml)
Expiration Date	12 months from date of receipt.



Immunohistochemistry staining of GPR88 i...