


Product Datasheet

GPR68 Antibody (orb1250491)

Description	GPR68 Antibody
Species/Host	Rabbit
Reactivity	Bat, Equine, Gorilla, Human, Primate
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	GPR68 antibody was raised against a peptide located in the N-Terminal domain of GPR68 (Human).
Target	GPR68
Preservatives	PBS, 0.1% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/ml
Storage	GPR68 antibody should be stored long term (months) at -80 °C and short term (days) at 4 °C. As with all antibodies avoid freeze/thaw cycles.
Note	For research use only
Application notes	GPR68 antibody can be used in ELISA, and immunohistochemistry starting at 10 µg/mL.
Clonality	Polyclonal
Uniprot ID	Q15743
NCBI	Q15743
Dilution Range	GPR68 antibody can be used in ELISA, and immunohistochemistry starting at 10 µg/mL.
Expiration Date	12 months from date of receipt.



Immunohistochemistry staining of GPR68 i...