


Product Datasheet

EPHA4 Antibody (orb1250514)

Description	EPHA4 Antibody
Species/Host	Rabbit
Reactivity	Bat, Bovine, Canine, Equine, Gorilla, Hamster, Human, Monkey, Mouse, Porcine, Primate, Rabbit, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	EPHA4 antibody was raised against a peptide located in the C-Terminal domain of EPHA4 (Human).
Target	EPHA4
Preservatives	PBS, 0.1% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/ml
Storage	EPHA4 antibody should be stored long term (months) at -80 °C and short term (days) at 4 °C. As with all antibodies avoid freeze/thaw cycles.
Note	For research use only
Application notes	EPHA4 antibody can be used in immunohistochemistry starting at 7.3 µg/mL.
Clonality	Polyclonal
Uniprot ID	P54764
NCBI	P54764
Dilution Range	EPHA4 antibody can be used in immunohistochemistry starting at 7.3 µg/mL.
Expiration Date	12 months from date of receipt.



Immunohistochemistry staining of EPHA4 i...