


## Product Datasheet

### SLC5A3 Antibody (orb1250516)

|                            |                                                                                                                                               |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | SLC5A3 Antibody                                                                                                                               |
| <b>Species/Host</b>        | Rabbit                                                                                                                                        |
| <b>Reactivity</b>          | Bat, Bovine, Canine, Hamster, Human, Monkey, Rabbit                                                                                           |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                  |
| <b>Tested Applications</b> | IHC                                                                                                                                           |
| <b>Immunogen</b>           | SLC5A3 antibody was raised against a peptide located in the C-Terminal domain of SLC5A3 (Human).                                              |
| <b>Target</b>              | SLC5A3                                                                                                                                        |
| <b>Preservatives</b>       | PBS, 0.1% sodium azide.                                                                                                                       |
| <b>Form/Appearance</b>     | Liquid                                                                                                                                        |
| <b>Concentration</b>       | 1 mg/ml                                                                                                                                       |
| <b>Storage</b>             | SLC5A3 antibody should be stored long term (months) at -80 °C and short term (days) at 4 °C. As with all antibodies avoid freeze/thaw cycles. |
| <b>Note</b>                | For research use only                                                                                                                         |
| <b>Application notes</b>   | SLC5A3 antibody can be used in ELISA, Western Blot, immunohistochemistry starting at 5 µg/mL, and immunofluorescence starting at 10 µg/mL.    |
| <b>Clonality</b>           | Polyclonal                                                                                                                                    |
| <b>Uniprot ID</b>          | <a href="#">P53794</a>                                                                                                                        |
| <b>NCBI</b>                | <a href="#">P53794</a>                                                                                                                        |
| <b>Dilution Range</b>      | SLC5A3 antibody can be used in ELISA, Western Blot, immunohistochemistry starting at 5 µg/mL, and immunofluorescence starting at 10 µg/mL.    |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                               |



Immunohistochemistry staining of SLC5A3 ...