

Product Datasheet

BCL2L15 antibody (orb125053)

Description	Goat polyclonal antibody to BCL2L15
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Target	BCL2L15
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Preliminary experiments in lysates of cell lines A549, HEK293 and Jurkat gave no specific signal but low background (at antibody concentration up to 1µg/ml).
Protein Sequence	SPTLQVASRN
Clonality	Polyclonal
MW	17.7
Entrez	440603
Dilution Range	ELISA: 1:8000, WB: 1 µg/ml), IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.