


Product Datasheet

BBC3 antibody (orb125054)

Description	Goat polyclonal antibody to BBC3
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	BBC3
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 24kDa band observed in Human peripheral blood mononucleocytes lysates (calculated MW of 26.5kDa according to NP_001120712.1). Recommended concentration: 0.1-0.3µg/ml.
Protein Sequence	DGGRPQRTVR
Clonality	Polyclonal
MW	26.5; 14.5; 10.3; 20.5
Entrez	27113
Dilution Range	ELISA: 1:16000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Peripheral Bloo...