


Product Datasheet

CD69 Antibody (FITC) (orb1250563)

Description	CD69 Antibody (FITC)
Species/Host	Hamster
Reactivity	Mouse
Conjugation	FITC
Tested Applications	FC, IHC, IP
Immunogen	Y245 mouse dendritic epidermal T cell line
Target	CD69
Preservatives	PBS, sodium azide.
Form/Appearance	Liquid
Concentration	0.5 mg/ml
Storage	Short term: +4°C; Long term: -20°C; Protect from light.
Note	For research use only
Application notes	CD69 antibody can be used in Western Blot starting at 0.5 - 4 µg/mL, immunohistochemistry starting at 2.5 µg/mL, and immunoprecipitation.
Isotype	IgG ₁
Clonality	Monoclonal
Uniprot ID	Q07108
NCBI	Q07108
Dilution Range	CD69 antibody can be used in Western Blot starting at 0.5 - 4 µg/mL, immunohistochemistry starting at 2.5 µg/mL, and immunoprecipitation.
Expiration Date	12 months from date of receipt.



Immunohistochemistry of mouse spleen tis...