

Product Datasheet

IFNG Antibody (PE) (orb1250575)

Description	IFNG Antibody (PE)
Species/Host	Mouse
Reactivity	Human, Monkey, Primate
Conjugation	PE
Tested Applications	FC
Immunogen	E. coli-expressed human IFN-gamma
Target	IFNG
Preservatives	Supplied in PBS/NaN ₃ and a stabilizing agent
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	IFN-gamma Antibody [B27] for use in flow cytometry assays.
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P01579
NCBI	P01579
Dilution Range	IFN-gamma Antibody [B27] for use in flow cytometry assays.
Expiration Date	12 months from date of receipt.