


Product Datasheet

HLA-DQA2 antibody (orb125058)

Description	Goat Polyclonal to HLA-DQA2. NP_064440.1; NP_002113.2, PMID: 20237496 and GeneID: 3118; 3117.
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	HLA-DQA2 & HLA-DQA1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 28kDa band observed in Human Bone Marrow lysates (calculated MW of 28.0kDa according to NP_064440.1). Recommended concentration: 1-3µg/ml.
Protein Sequence	QGLRSVGASRH
Clonality	Polyclonal
MW	28
Entrez	3117 , 3118
Dilution Range	ELISA: 1:32000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Bone Marr...