

Product Datasheet

FAS Antibody (Biotin) (orb1250603)

Description	FAS Antibody (Biotin)
Species/Host	Mouse
Reactivity	Human, Monkey, Primate
Conjugation	Biotin
Tested Applications	Activity Assay, Blocking, COSTIM, FC, IP, Stim
Immunogen	Human CD95 transfected L cells
Target	FAS
Preservatives	Supplied in PBS/NaN ₃
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD95 Antibody [DX2] for use in flow cytometry, immunohistochemistry / immunocytochemistry, immunoprecipitation, and ELISA assays.
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P25445
NCBI	P25445
Dilution Range	CD95 Antibody [DX2] for use in flow cytometry, immunohistochemistry / immunocytochemistry, immunoprecipitation, and ELISA assays.
Expiration Date	12 months from date of receipt.