


Product Datasheet

CD79B Antibody (APC) (orb1250606)

Description	CD79B Antibody (APC)
Species/Host	Mouse
Reactivity	Human
Conjugation	APC
Tested Applications	FC
Immunogen	Ramos B cells
Target	CD79B
Preservatives	Supplied in PBS/NaN3 and a stabilizing agent
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD79b Antibody [CB3-1] for use in flow cytometry assays.
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P40259
NCBI	P40259
Dilution Range	CD79b Antibody [CB3-1] for use in flow cytometry assays.
Expiration Date	12 months from date of receipt.



Human
peripheral
blood
lymphocytes
were ...