

Product Datasheet

CD55 Antibody (PE/Cy5) (orb1250620)

Description	CD55 Antibody (PE/Cy5)
Species/Host	Mouse
Reactivity	Human
Conjugation	PE/Cy5
Tested Applications	Blocking, FC, IHC-Fr, IP
Immunogen	PHA activated peripheral blood mononuclear cells
Target	CD55
Preservatives	Supplied in PBS/NaN ₃ and a stabilizing agent
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD55 Antibody [143-30] for use in flow cytometry assays.
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P08174
NCBI	P08174
Dilution Range	CD55 Antibody [143-30] for use in flow cytometry assays.
Expiration Date	12 months from date of receipt.