

Product Datasheet

CD55 Antibody (FITC) (orb1250622)

Description	CD55 Antibody (FITC)
Species/Host	Mouse
Reactivity	Human
Conjugation	FITC
Tested Applications	FC
Immunogen	PHA activated peripheral blood mononuclear cells
Target	CD55
Preservatives	Supplied in PBS/NaN ₃
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD55 Antibody [143-30] for use in flow cytometry, immunohistochemistry, and western blotting assays.
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P08174
NCBI	P08174
Dilution Range	CD55 Antibody [143-30] for use in flow cytometry, immunohistochemistry, and western blotting assays.
Expiration Date	12 months from date of receipt.