

Product Datasheet

PTPRC Antibody (APC) (orb1250634)

Description	PTPRC Antibody (APC)
Species/Host	Mouse
Reactivity	Equine, Human
Conjugation	APC
Tested Applications	FC
Immunogen	Purified T cells from human lymph nodes
Target	PTPRC
Preservatives	Supplied in PBS/NaN ₃ and a stabilizing agent
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD45 Antibody [F10-89-4] for use in flow cytometry assays.
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	P08575
NCBI	P08575
Dilution Range	CD45 Antibody [F10-89-4] for use in flow cytometry assays.
Expiration Date	12 months from date of receipt.