


Product Datasheet

PTPRC Antibody (PE/Texas Red) (orb1250635)

Description	PTPRC Antibody (PE/Texas Red)
Species/Host	Mouse
Reactivity	Equine, Human
Conjugation	PE/Texas Red
Tested Applications	FC
Immunogen	Purified T cells from human lymph nodes
Target	PTPRC
Preservatives	Supplied in PBS/NaN3 and a stabilizing agent
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD45 Antibody [F10-89-4] for use in flow cytometry assays.
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	P08575
NCBI	P08575
Dilution Range	CD45 Antibody [F10-89-4] for use in flow cytometry assays.
Expiration Date	12 months from date of receipt.



Human peripheral blood lymphocytes were ...