

Product Datasheet

MPV17 antibody (orb125064)

Description

Goat polyclonal antibody to MPV17

Species/Host

Goat

Reactivity

Bovine, Canine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications
Target

MPV17

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

WB: Approx 19kDa band observed in Human Brain (Frontal Cortex) lysates (calculated MW of 19.7kDa according to NP_002428.1). Recommended concentration: 1-3µg/ml. An additional band of unknown identity was also consistently observed at 26kDa. This band was succ. IHC: In paraffin embedded Human Liver shows heavy textured cytoplasm staining in hepatocytes, consistent with mitochondria.. Recommended concentration, 5-10µg/ml.

Protein Sequence

QRALAAHPWKVQ

Clonality

Polyclonal

MW


19.7

Entrez
[4358](#)
Dilution Range


ELISA: 1:32000, WB: 1-3 µg/ml, IHC-P: 5-10µg/ml

Expiration Date

12 months from date of receipt.

Western blot analysis
of Human Frontal C...



Immunohistochemical
staining of Human
Li...