


Product Datasheet

CD33 Antibody (FITC) (orb1250658)

Description	CD33 Antibody (FITC)
Species/Host	Mouse
Reactivity	Human
Conjugation	FITC
Tested Applications	FC
Target	CD33
Preservatives	Supplied in PBS/NaN3
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD33 Antibody [WM53] for use in flow cytometry, immunohistochemistry / immunocytochemistry, and western blotting assays.
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P20138
NCBI	P20138
Dilution Range	CD33 Antibody [WM53] for use in flow cytometry, immunohistochemistry / immunocytochemistry, and western blotting assays.
Expiration Date	12 months from date of receipt.



Human peripheral blood lymphocytes were ...