



Product Datasheet

ALDH6A1 antibody (orb125068)

Description	Goat polyclonal antibody to ALDH6A1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	ALDH6A1 (aa487-496)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 55kDa band observed in Human, Mouse and Rat Kidney lysates (calculated MW of 57.8kDa according to NP_005580.1). Recommended concentration: 0.1-0.3µg/ml.. IHC: In paraffin embedded Human Liver shows textured staining of the cytoplasm consistent with mitochondria. Recommended concentration, 5-10µg/ml.
Protein Sequence	SRSSFRGDTN
Clonality	Polyclonal
MW	57.8; 56.2; 41.2
Entrez	4329
Dilution Range	ELISA: 1:32000, WB: 0.1-0.3 µg/ml, IHC-P: 5-10µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse Kidney ly...



Immunohistochemical staining of Human Li...