

Product Datasheet

TRIM72 antibody (orb125071)

Description

Goat Polyclonal to TRIM72. NP_001008275.1, PMID:

Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	TRIM72
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 52kDa band observed in Mouse Heart lysates (calculated MW of 52.8kDa according to Mouse NP_001073401.1). Recommended concentration: 0.01-0.03µg/ml. Primary incubation was 1 hour.. IHC: In paraffin embedded Human Skeletal Muscle shows staining of the fibre membranes, detached artificially by the sectioning. Recommended concentration: 3-5µg/ml.
Protein Sequence	RPQALSTNLQ
Clonality	Polyclonal
MW	52.7
Entrez	493829
Dilution Range	ELISA: 1:16000, WB: 0.01-0.03 µg/ml, IHC-P: 3-5 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human Skeletal Muscle



Western blot analysis of Mouse Heart lysates