

Product Datasheet

CD4 Antibody (PE/Texas Red) (orb1250719)

Description	CD4 Antibody (PE/Texas Red)
Species/Host	Mouse
Reactivity	Human
Conjugation	PE/Texas Red
Tested Applications	FC
Immunogen	Thymocytes and E-rosetted lymphocytes
Target	CD4
Preservatives	Supplied in PBS/NaN ₃ and a stabilizing agent
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD4 Antibody [RFT4] for use in flow cytometry assays.
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P01730
NCBI	P01730
Dilution Range	CD4 Antibody [RFT4] for use in flow cytometry assays.
Expiration Date	12 months from date of receipt.