


Product Datasheet

ITGA4 Antibody (PE) (orb1250751)

Description	ITGA4 Antibody (PE)
Species/Host	Mouse
Reactivity	Human
Conjugation	PE
Tested Applications	FC, IHC-Fr, IHC-P
Target	ITGA4
Preservatives	Supplied in PBS/NaN3 and a stabilizing agent
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD49d Antibody [7.2] for use in flow cytometry assays.
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P13612
NCBI	P13612
Dilution Range	CD49d Antibody [7.2] for use in flow cytometry assays.
Expiration Date	12 months from date of receipt.



Human
peripheral
blood
lymphocytes
were ...