

## Product Datasheet

### ITGA4 Antibody (FITC) (orb1250752)

|                            |                                                                                                   |
|----------------------------|---------------------------------------------------------------------------------------------------|
| <b>Description</b>         | ITGA4 Antibody (FITC)                                                                             |
| <b>Species/Host</b>        | Mouse                                                                                             |
| <b>Reactivity</b>          | Human                                                                                             |
| <b>Conjugation</b>         | FITC                                                                                              |
| <b>Tested Applications</b> | FC                                                                                                |
| <b>Target</b>              | ITGA4                                                                                             |
| <b>Preservatives</b>       | Supplied in PBS/NaN <sub>3</sub>                                                                  |
| <b>Form/Appearance</b>     | Liquid                                                                                            |
| <b>Concentration</b>       | batch dependent                                                                                   |
| <b>Storage</b>             | Store vial at 2-8°C                                                                               |
| <b>Note</b>                | For research use only                                                                             |
| <b>Application notes</b>   | CD49d Antibody [7.2] for use in flow cytometry, immunocytochemistry, and western blotting assays. |
| <b>Isotype</b>             | IgG1                                                                                              |
| <b>Clonality</b>           | Monoclonal                                                                                        |
| <b>Uniprot ID</b>          | <b>P13612</b>                                                                                     |
| <b>NCBI</b>                | <b>P13612</b>                                                                                     |
| <b>Dilution Range</b>      | CD49d Antibody [7.2] for use in flow cytometry, immunocytochemistry, and western blotting assays. |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                   |