


Product Datasheet

PTCHD1 antibody (orb125076)

Description	Goat Polyclonal to PTCHD1. NP_775766.2, PMID: 20531469 and GenelD: 139411.
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	PTCHD1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 100kDa band observed in Human Brain (Cerebellum) lysates (calculated MW of 101kDa according to NP_775766.2). Recommended concentration: 0.5-2µg/ml. An additional band of unknown identity was also consistently observed at 50kDa. This band was succes
Protein Sequence	DFQKTSRVSER
Clonality	Polyclonal
MW	101
Entrez	139411
Dilution Range	ELISA: 1:16000, WB: 0.5-2 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Cerebellu...