


Product Datasheet

RBM20 antibody (orb125077)

Description	Goat Polyclonal to RBM20. NP_001127835.1, PMID: 20590677 and GenelD: 282996.
Species/Host	Goat
Reactivity	Bovine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	RBM20 (aa717-727)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 150kDa band observed in lysates of cell line HEK293 (calculated MW of 134kDa according to Human NP_001127835.2). Preliminary experiments also showed a 150kDa band in IU2OS cell lysate and in Human Smooth Muscle lysate. This molecular weight is rout
Protein Sequence	QLDKAELDERP
Clonality	Polyclonal
MW	134
Entrez	282996
Dilution Range	ELISA: 1:64000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse Heart lys...