


## Product Datasheet

### CD19 Antibody (PE/Texas Red) (orb1250779)

|                            |                                                           |
|----------------------------|-----------------------------------------------------------|
| <b>Description</b>         | CD19 Antibody (PE/Texas Red)                              |
| <b>Species/Host</b>        | Mouse                                                     |
| <b>Reactivity</b>          | Human                                                     |
| <b>Conjugation</b>         | PE/Texas Red                                              |
| <b>Tested Applications</b> | FC                                                        |
| <b>Immunogen</b>           | NALM-1 and NALM-16 leukemia cell line                     |
| <b>Target</b>              | CD19                                                      |
| <b>Preservatives</b>       | Supplied in PBS/NaN3 and a stabilizing agent              |
| <b>Form/Appearance</b>     | Liquid                                                    |
| <b>Concentration</b>       | batch dependent                                           |
| <b>Storage</b>             | Store vial at 2-8°C                                       |
| <b>Note</b>                | For research use only                                     |
| <b>Application notes</b>   | CD19 Antibody [SJ25-C1] for use in flow cytometry assays. |
| <b>Isotype</b>             | IgG1                                                      |
| <b>Clonality</b>           | Monoclonal                                                |
| <b>Uniprot ID</b>          | <b>P15391</b>                                             |
| <b>NCBI</b>                | <b>P15391</b>                                             |
| <b>Dilution Range</b>      | CD19 Antibody [SJ25-C1] for use in flow cytometry assays. |
| <b>Expiration Date</b>     | 12 months from date of receipt.                           |



Human  
peripheral  
blood  
lymphocytes  
were ...