


Product Datasheet

ITGAL Antibody (Biotin) (orb1250786)

Description	ITGAL Antibody (Biotin)
Species/Host	Mouse
Reactivity	Human
Conjugation	Biotin
Tested Applications	FC
Immunogen	Fibronectin purified monocytes
Target	ITGAL
Preservatives	Supplied in PBS/NaN3
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD11a Antibody [38] for use in flow cytometry, immunohistochemistry, ELISA, and immunoprecipitation assays.
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	P20701
NCBI	P20701
Dilution Range	CD11a Antibody [38] for use in flow cytometry, immunohistochemistry, ELISA, and immunoprecipitation assays.
Expiration Date	12 months from date of receipt.



Human peripheral blood lymphocytes were ...