

Product Datasheet

TNFSF10 Antibody (orb1250836)

Description	TNFSF10 Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-Fr, IP, WB
Immunogen	Recombinant TRAIL
Target	TNFSF10
Preservatives	Supplied in borate buffered saline, pH 8.2; No preservatives or amine-containing buffer salts added
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	TRAIL Antibody [SB91c] for use in ELISA assays.
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	P50591
NCBI	P50591
Dilution Range	TRAIL Antibody [SB91c] for use in ELISA assays.
Expiration Date	12 months from date of receipt.