

Product Datasheet

TNFSF10 Antibody (Biotin) (orb1250837)

Description	TNFSF10 Antibody (Biotin)
Species/Host	Mouse
Reactivity	Human
Conjugation	Biotin
Tested Applications	FC
Immunogen	Recombinant TRAIL
Target	TNFSF10
Preservatives	Supplied in PBS/NaN ₃
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	TRAIL Antibody [SB91a] for use in ELISA assays.
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	P50591
NCBI	P50591
Dilution Range	TRAIL Antibody [SB91a] for use in ELISA assays.
Expiration Date	12 months from date of receipt.