



Product Datasheet

Influenza A, Nucleoprotein Antibody [FluA-NP 4F1] (orb1250841)

Description	Influenza A, Nucleoprotein Antibody [FluA-NP 4F1]
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Recombinant influenza virus type A nucleoprotein
Preservatives	Supplied in borate buffered saline, pH 8.2; No preservatives or amine-containing buffer salts added
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	Influenza A, Nucleoprotein Antibody [FluA-NP 4F1] for use in immunocytochemistry, ELISA, and western blotting assays.
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	Influenza A, Nucleoprotein Antibody [FluA-NP 4F1] for use in immunocytochemistry, ELISA, and western blotting assays.
Expiration Date	12 months from date of receipt.



Recombinant influenza proteins were reso...



MDCK cells infected with the influenza A...