


Product Datasheet

Csf2 Antibody (Biotin) (orb1250846)

Description	Csf2 Antibody (Biotin)
Species/Host	Rat
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	FC
Immunogen	COS-expressed mouse GM-CSF
Target	Csf2
Preservatives	Supplied in PBS/NaN3
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	GM-CSF Antibody [MP1-31G6] for use in immunohistochemistry / immunocytochemistry, ELISA, immunoprecipitation, and multiplex assays.
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P01587
NCBI	P01587
Dilution Range	GM-CSF Antibody [MP1-31G6] for use in immunohistochemistry / immunocytochemistry, ELISA, immunoprecipitation, and multiplex assays.
Expiration Date	12 months from date of receipt.



Standard curve generated with Rat Anti-M...