



Product Datasheet

CYP24A1 antibody (orb125085)

Description	Goat Polyclonal to CYP24A1. NP_000773.2; NP_001122387.1, PMID: 20847308 and GenElD: 1591.
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	CYP24A1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 50kDa band observed in Human and Rat Liver lysates (calculated MW of 51.2kDa according to NP_001122387.1). Recommended concentration: 0.1-0.3µg/ml. Primary incubation was 1 hour.. IHC: Paraffin embedded Human Skin. Recommended concentration: 3.75µg/ml.
Protein Sequence	RKYDIQATDN
Clonality	Polyclonal
MW	58.9; 51.2
Entrez	1591
Dilution Range	ELISA: 1:32000, WB: 0.1-0.3 µg/ml, IHC-P: 3.75 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Liver lys...



Immunohistochemical staining of Human Sk...