


Product Datasheet

IL5 Antibody (Biotin) (orb1250874)

Description	IL5 Antibody (Biotin)
Species/Host	Rat
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	Blocking, Depletion, FC, ICC, IHC-Fr, IP
Immunogen	Semi-purified T cell clone supernatant
Target	IL5
Preservatives	Supplied in PBS/NaN3
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	IL-5 Antibody [TRFK4] for use in immunohistochemistry / immunocytochemistry, ELISA, ELISpot, and multiplex assays.
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	P04401
NCBI	P04401
Dilution Range	IL-5 Antibody [TRFK4] for use in immunohistochemistry / immunocytochemistry, ELISA, ELISpot, and multiplex assays.
Expiration Date	12 months from date of receipt.



Standard curve generated with Rat Anti-M...