


Product Datasheet

PCBP1 antibody (orb125088)

Description	Goat Polyclonal to PCBP1. CAA55016.1, PMID: 19211566 and GeneID: 5093.
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	PCBP1 (aa234-247)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 38kDa band observed in lysates of cell line HeLa and in Human Placenta lysates (calculated MW of 37.5kDa according to NP_006187.2). Recommended concentration: 0.01-0.03µg/ml.. IHC: Paraffin embedded Human Lung (respiratory epithelium). Recommended concentration: 5µg/ml.
Protein Sequence	ADLAKLNQVARQQSH
Clonality	Polyclonal
MW	37.5
Entrez	5093
Dilution Range	ELISA: 1:32000, WB: 0.01-0.03 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa lysate using PCBP1 antibody (orb125088).